

Product datasheet for **TP504614**

Fam166a (NM_026624) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse family with sequence similarity 166, member A (Fam166a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR204614 representing NM_026624 Red =Cloning site Green =Tags(s) MTATQKHNLFTPEPHYIPGYAGFYPLRLYQVGNTYGRTTAQLLTDPSPVQKSPCSVLSPMTKPKFIEDFSK SKPPWIPCRDLREPYIPHYTSLKPYKNFEILGQLPRQDQDVTQGPQVENRQGPLTAGFMPYPYPACPPG RKGEARDLGHPGLLLAYGEEAWKDAAPLQDTPGKNNQLYHCRDEYLPPHPPQETLDVGRFQRLPQLD HP NLIQRKAISGYAGFVPRFAWVMGMNYRDGVTQAMDDFDKNQFVFRHPVCALGERLPRTWHPNTTIYRS QG LIPFYMGMFIPSMQDNALTFGNSTRRAYQKELDRRSHTL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	36.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_080900</u>
Locus ID:	68222


[View online »](#)

UniProt ID:	<u>Q9D4K5</u>
RefSeq Size:	1130
Cytogenetics:	2 A3
RefSeq ORF:	957
Synonyms:	4931415M17