

Product datasheet for **TP504650**

Lgals9 (NM_001159301) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse lectin, galactose binding, soluble 9 (Lgals9), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA	>MR204650 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MALFSAQSPYINPIIPFTGPIQGGGLQEGLQVTLQGTTKSFAQRFVWNFQNSFNGNDIAFHFNPRFEEGGY
VVCNTKQNGQWGPEERKMQMPFQKGMPELFCFLVQRSEFKVMVNKKFFVQYQHRVPYHLVDTIAVSGCLK
LSFITFQTQDFRPAHQAPMAQTTHMVHSTPGQMFSTPGIPPVYPTPAYTIPFYTPIPNGLYPSKSIMI
SGNVLPDATRFHINLRCGGDIAFHNLNPRFNENAVVRNTQINNSWGQEERSLLGRMPFSRGQSFSVWIICE
GHCFKVAVNGQHMCEYYHRLKNLQDINTLEVAGDIQLTHVQT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	36.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001152773
Locus ID:	16859
UniProt ID:	O08573



[View online »](#)

RefSeq Size:	1493
Cytogenetics:	11 B5
RefSeq ORF:	969
Synonyms:	AA407335; AI194909; AI265545; gal-9; galectin-9; Lgals5; LGALS35
Summary:	<p>Binds galactosides (By similarity). Has high affinity for the Forssman pentasaccharide (By similarity). Ligand for HAVCR2/TIM3 (By similarity). Binding to HAVCR2 induces T-helper type 1 lymphocyte (Th1) death (By similarity). Also stimulates bactericidal activity in infected macrophages by causing macrophage activation and IL1B secretion which restricts intracellular bacterial growth (PubMed:20937702). Ligand for P4HB; the interaction retains P4HB at the cell surface of Th2 T-helper cells, increasing disulfide reductase activity at the plasma membrane, altering the plasma membrane redox state and enhancing cell migration (PubMed:21670307). Ligand for CD44; the interaction enhances binding of SMAD3 to the FOXP3 promoter, leading to up-regulation of FOXP3 expression and increased induced regulatory T (iTreg) cell stability and suppressive function (PubMed:25065622). Promotes ability of mesenchymal stromal cells to suppress T-cell proliferation (By similarity). Expands regulatory T-cells and induces cytotoxic T-cell apoptosis following virus infection (By similarity). Activates ERK1/2 phosphorylation inducing cytokine (IL-6, IL-8, IL-12) and chemokine (CCL2) production in mast and dendritic cells (By similarity). Inhibits degranulation and induces apoptosis of mast cells (By similarity). Induces maturation and migration of dendritic cells (By similarity). Inhibits natural killer (NK) cell function (PubMed:23408620). Can transform NK cell phenotype from peripheral to decidual during pregnancy (By similarity). Astrocyte derived galectin-9 enhances microglial TNF production (PubMed:25158758). May play a role in thymocyte-epithelial interactions relevant to the biology of the thymus. May provide the molecular basis for urate flux across cell membranes, allowing urate that is formed during purine metabolism to efflux from cells and serving as an electrogenic transporter that plays an important role in renal and gastrointestinal urate excretion (By similarity). Highly selective to the anion urate (By similarity).[UniProtKB/Swiss-Prot Function]</p>