

Product datasheet for **TP504804**

Bcl2l14 (NM_025778) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse BCL2-like 14 (apoptosis facilitator) (Bcl2l14), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR204804 protein sequence Red =Cloning site Green =Tags(s)
	<p>MCSTSVYDLEDIPLDDDDPNSIEFKILAFYARHHVFNTPAVFSPKLSRTRLSQKALGTWSTDSWTQVS LPCRGSPSEKNISLGKKKSSWRTLFRVAEKEEGPPSPKEIRAQGPQGPFVVERQSGFHNQHWPRSLSS VEQRLESEVWDSKVACIANRVAEIVSWPPDPVIHSQGGSKLKERVSEILYFRFEGPCDSKKDGEDQIIS KIVELLKFSGDQLGREIKKDKALMSSFQDGLSYSTFKTITDLFLRDVDRGESEVKAQGFKAALIDAIA KLTAIDNHPMNRMLGFGTKYLREYFSPWVQQNGGWEEKILGISHEEVD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	36.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080054
Locus ID:	66813
UniProt ID:	Q9CPT0



[View online »](#)

RefSeq Size: 2220
Cytogenetics: 6 G1
RefSeq ORF: 984
Synonyms: 4930452K23Rik; 4933405K19Rik; 9030625M01Rik; AI429211; Bclg
Summary: Plays a role in apoptosis.[UniProtKB/Swiss-Prot Function]