

## Product datasheet for TP505012

### Snx15 (NM\_026912) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse sorting nexin 15 (Snx15), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205012 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MSRQAKDDFLRHYTVSDPRTHPKGYTEYKVTAQFISKKDPEDIKEVWWKRYSDFRKLHGDLAYTHRNLF RRLEEFPAFPRAQVFGREFASVIEERRKGAEDLLRFTVPIPALNNSPQLKEFFRGGEVTRPSEVSRDLRI LPPPLIPTPPPDEARLLQPLPAERRGQEELEVPVDPLPSSPAQEALDLLFSCDSTEEASSSLARGPLSEA ELALFDPYSKEESTGPSPTHTGELAAIEVESKRLDQEPWEPGGQEEEEAEDGEPAPAYLGQATELITQAL RNEKAGAYAAALQGYQDGVHILLQGVSGDPSPARREGVKKKAAEYLKRAEMLHHTLP  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-MYC/DDK
Predicted MW:	37.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_081188</a></u>
Locus ID:	69024
UniProt ID:	<u><a href="#">Q91WE1</a></u>


[View online »](#)

RefSeq Size:	1554
Cytogenetics:	19 A
RefSeq ORF:	1011
Synonyms:	1500032B08Rik; E130013C21Rik
Summary:	May be involved in several stages of intracellular trafficking. Overexpression of SNX15 disrupts the normal trafficking of proteins from the plasma membrane to recycling endosomes or the TGN (By similarity).[UniProtKB/Swiss-Prot Function]