

Product datasheet for **TP505055**

Anxa9 (NM_001085383) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse annexin A9 (Anxa9), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205055 protein sequence Red =Cloning site Green =Tags(s)
	MSASCGPLGTSLTQEILSSLGLADKTAAWGTLGTLRTFLSFSVDKDVQRLLKAIAGQGVDYDTIVDVLTN RSREQRQLISRAFQERTKQDLLKSLQAALSGNLEKIVVALLQPAAQFDAQELRTALKTSGSAEDVALEIL ATRAAPGLQACLA VYKHDFQVEAEEDIRTETNGILQDLLLALSKGDRESYSGIIDYNLEEQDVRLQQAG ESSTAGQWVLLLTQRSPEHLIRVFDQYRRCTGQELEDAIRNCFHGDAQLALISLASMLRNTALYFANKLH QALQETEPNFQVLTRVLISRSEDLISIRAEFKKFGKSLYSSLQDVVRGDCRSALLALCRAEDI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	38.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



RefSeq:	<u>NP_001078852</u>
Locus ID:	71790
UniProt ID:	<u>Q9JHQ0</u>
RefSeq Size:	1895
Cytogenetics:	3 F2.1
RefSeq ORF:	1035
Synonyms:	1110003P15Rik; 2310069F17Rik
Summary:	May act as a low affinity receptor for acetylcholine.[UniProtKB/Swiss-Prot Function]