

Product datasheet for **TP505101**

Armcl2 (NM_026290) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse armadillo repeat containing 12 (Armcl2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205101 representing NM_026290 Red=Cloning site Green=Tags(s) MGKTIPRFLEQLDIKSFVGLATGAGALYLLYKAVRTGLKCHPPLCSNSPICIARLAIERERHGRDSGEI RRLNSLDCKQDEYTRSMILHNITRCVYLLEAEASSCTMDDIDLVDMLDEKDNSVKIQALNALKAFSGI RKFRLLKIQEHCIKVLLEISTIWDLELHVAGLRLLNNLPLPDYVHPQLRRVMPALMEIIQSDCILAQVQAV RLLSYLAQKNDLLYDILNCQVRPNFLNFQSSQPGSLLFEVLVFAEHLSEGRNATHYRAVKWHYNEQSLH EALFGDESRLADRLLSLVHPPEEVQIQACKVIVSLQCPQDLGSRPSSCRPSHSCFKTGK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	38.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



RefSeq:	<u>NP_080566</u>
Locus ID:	67645
UniProt ID:	<u>Q80X86</u>
RefSeq Size:	1103
Cytogenetics:	17 A3.3
RefSeq ORF:	1020
Synonyms:	4930511111Rik; AV041658