

Product datasheet for TP505203

Ahsg (NM_013465) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse alpha-2-HS-glycoprotein (Ahsg), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Tag:	C-MYC/DDK
Predicted MW:	37.3 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_038493
Locus ID:	11625
RefSeq Size:	1524
Cytogenetics:	16 13.74 cM
RefSeq ORF:	1038
Synonyms:	fetuin-A
Summary:	Probably involved in differentiation. [UniProtKB/Swiss-Prot Function]



[View online »](#)