

Product datasheet for TP505218

Lrp2bp (NM_026278) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse Lrp2 binding protein (Lrp2bp), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205218 representing NM_026278 Red =Cloning site Green =Tags(s) MKLSSEKLPKNPFSLSQYAAKQQKFFQWKKEKPDYYLHANLVDTALQFLKERIRRGDAMAYFLRGQLYFE EGWYEEALAQFEEIQEKDHQAIYQLGVMMYDGLGTIANAEKGVNMRKILDSSCPQTMHLKFAAAYNLGR AYFEGKGVKRSDEEAERLWLLAADNGNPKASVKAQSILGLFYSMKPEKELEKAFFWHSEACGNGSLESQG ALGLMYFYGGQIRQDTDAALHCLREAAERGNVYAQGTLVEYYYYKMKFFTCKVSFSKRIADYDEVHDIPMI AHVTDCLPEFIIKGMAMAAFYHGRCLQLGLGIMKDEESAKHYYSKACRLNPTLADELHSLIHQRI TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	40.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080554
Locus ID:	67620
UniProt ID:	Q9D4C6


[View online »](#)

RefSeq Size:	3420
Cytogenetics:	8 B1.1
RefSeq ORF:	1038
Synonyms:	1700113N17Rik; 2310006J04Rik; 4930479L12Rik; MegBP
Summary:	May act as an adapter that regulates LRP2 function.[UniProtKB/Swiss-Prot Function]