

Product datasheet for **TP505335**

Fam118b (NM_194257) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse family with sequence similarity 118, member B (Fam118b), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205335 protein sequence Red =Cloning site Green =Tags(s) MASTGSQASDKNKIFGFFNDGEPPTKKPRKLLPSLKTKKPRELVLVIGTGISAAPQVPALKSWKGLIQ ALLDAAIDFDLLEDEESKKFQKCLHEDKNLVHVAHDLIQKLSPTSNVRSTFFKDCLYEVFDDLESKMED SGKQLLQSVLHLMENGALVLTNFDNLLLEYAADQKGKQLESGLTDEKKVLEWAQEKRLSVLHIHGVYT NPSGIVLHPAGYQNVLRNTEVMREIQKLYENKSFLFLGCGWTVDDTTFQALFLEAVKHKSDLEHFMLVRR GDVDEFKKLRENMLDKGIKVISYGNDYADLPEYFKRLTCEISTRGRSAGMAREGQLNGSSAAHGEIRGCS T TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	39.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_919233
Locus ID:	109229
UniProt ID:	Q8C569


[View online »](#)

RefSeq Size:	1915
Cytogenetics:	9 A4
RefSeq ORF:	1053
Synonyms:	2310022O21Rik; 2700018L24Rik; C030004A17Rik; Ppp4r1l
Summary:	May play a role in Cajal bodies formation.[UniProtKB/Swiss-Prot Function]