

## Product datasheet for **TP505577**

### Ssbp4 (NM\_133772) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse single stranded DNA binding protein 4 (Ssbp4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205577 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MYGKASKGCAPSDGQAREKLALVYVEYLLHVGAQKSAQTFNSEIRWEKNITLGEPPGFLHSWWCVFWDLY CAAPDRREACEHSSEAKVFQDYSAAAAPSPVMGAMTPNDAMAAGPVAPGFFQPFMSPRFPGGARPILRMP GQPPVGLPGSQLIPAAMDPSRVQGHPSLGGPMQRVTPPRGMASVGPQGYGTGMRPPNSLATSQVLPS MNMGPGVGRGPWASPSGNSIPYSSSSPGSYSGPAGVGGAPGTPIMPSPGDSTNSENMYTIMNPIGPGAGR ANFPLGPSPEGPMASMSTMEPHHVNGSLGSGDMDGLPKNSPGAVGGLSNAPGTPRDDGEMAAAGTFLHPF PSESYPGMTMSV  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	37.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_598533</a>
Locus ID:	76900
UniProt ID:	<a href="#">Q3U4B1</a>



[View online »](#)

RefSeq Size: 1498  
Cytogenetics: 8 B3.3  
RefSeq ORF: 1092  
Synonyms: 1210002E11Rik; AW743380; Sspb4