

Product datasheet for TP505971

Klhdc3 (NM_001163729) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse kelch domain containing 3 (Klhdc3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205971 protein sequence Red=Cloning site Green=Tags(s) MLRWTVHLEGGPRRVNHA AVAVGHRVYSFGGYCSGEDYETLRQIDVHIFNAVSLRWTKLPPVRPAVRGQ A PVVPYMRYPGHSTVLIDDTVFLWGGRNDTEGACNVLYAFDVNTHKWSTPRVSGTVPGARDGHSACVLGKI M YIFGGYEQLADCFNSNDIHKLDTSTMTWTLVCTKGNPARWRDFHSATMLGNHMYVFGGRADRFPGPFHSN NE IYCNRIRVFDTRTEAWLDCPHTPVLPEGRRSHSAFGYNGELYIFGGYNARLNRHFHDLWKFNPGSFTWKK IEPKGKGPCPRRRQCCCIVGDKIVLFGGTSPSPEEGLGDEFDLIDHSDLHILDFSPSLKTLCKLAVIQYS LDQSCLPHDIRWELNAMTTNSNISRPVSSHG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	43.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: NP_001157201

Locus ID: 71765

UniProt ID: Q8VEM9

RefSeq Size: 1973

Cytogenetics: 17 C

RefSeq ORF: 1146

Synonyms: 1300011D16Rik; Peas

Summary: May be involved in meiotic recombination process.[UniProtKB/Swiss-Prot Function]