

Product datasheet for **TP505986**

Col9a3 (BC030945) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse procollagen, type IX, alpha 3 (cDNA clone MGC:32160 IMAGE:5003028), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR205986 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPGKDGKDGVPGLDGEKEAGRNGGQGEKGPNGLPGLPGRAGSKGEKGEPGRTGELGEAGPSGEPGIPGD</p> <p>VGVPGERGEAGHRGSVGALGPQGPPGAPGIRGFQGGKSTGDPGLPGPQGLRGDVGDRPGGATGPKGDQ</p> <p>GIAGSDGLPGDKGELGPNGPVGQKGESGSRGELGPKGIQGPNGTSGVQGVPGPPGLGLQGVQGVPGITG</p> <p>KPGVPGKEASEQRIRELCGGMISEQIAQLAAHLRKPLAPGSIGRPGPAGPPGPPGPPGSIGHPGARGPPGYRGPTGELGDPGPRGSQGDRGDKGATGAGLDGPAGDQGYQGPQGVPGISKDGRDGAHGEPGLPGDPGLPG</p> <p>AAGAQTGTPGICDTSACQGAVLGGGGEKSGPRSS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	36 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	12841
RefSeq Size:	1605
Cytogenetics:	2 103.14 cM
RefSeq ORF:	1149
Synonyms:	AV006866