

Product datasheet for **TP506064**

Cuedc1 (NM_198013) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse CUE domain containing 1 (Cuedc1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206064 protein sequence Red =Cloning site Green =Tags(s)
	MTSLFRRSSSGSGGGGATGARGAGTGAGDGSTAPQELNNSRPARQVRRLEFNQAMDDFKTMFPNMDY DII ECVLRANS GAVDATIDQLLQMNLEAGGGSAYEDSSDSIPPEILERTLEPDSSEEEPPPVYSPPAYHM HVFDRPYLMAPPTPPPRIDVPGSGQPASQRRYRNWNPPLLGSLPDDFLRILPQQMDSIQGHPGGSKPM SG EGGPPAPGPMACDQDSRWKQYLEDERIALFLQNEEFMKELQRNRDFFLLALERDLKYESQKSKSNYAAV GNDGGFPSSVPGTSETNPTVSEDALFRDKLKHMGKSTRRLKFELARAFSEKTKMRKSKKKHLPKLQSLGA AASTANLLDDVEGHAYEEDFRGRRQEVPKVEEALREGQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	42.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_932130</u>



[View online »](#)

Locus ID: 103841

UniProt ID: [Q8R3V6](#)

RefSeq Size: 3189

Cytogenetics: 11 C

RefSeq ORF: 1164

Synonyms: AI841487; C330016O16Rik