

Product datasheet for TP506234

Serpinb10 (NM_001160307) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse serine (or cysteine) peptidase inhibitor, clade B (ovalbumin), member 10 (Serpinb10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206234 protein sequence Red =Cloning site Green =Tags(s) MASLAVSINQFALEFSKKLAESAEGRNIFSPWGISTALAMLYLGTKGTTADQMAQVLQFSSVEDFKSCP DSEKKRKMEFNSGKFEEIQSDFQTLAAEILKPGNSYVLKTANRIYGEKTYPFHNKYLEDMKTYFGAEPQS VNFVEASGQIRKEINSWVGSQTGGKIPNLLPDDSVDTKTKMVLVNALYFKGTWEHQFSVKNTTERPFRVN KTTSKPVQMMSMKQSLQVFHIEELQTIGLQLHYQNRDLSLLLLPEAIDGLEQLERAITYEKLDKWTSAD MMDTYEVQLYLPKFKMEESYDLKSALRGMGMTDVSQSKADFSNMTSERNFLSNVFKHTFLEINEEGTE AAAGTGSEISVRIKAPSIELNVDHPFLFFIRHNKTKSILFCGRFCSP TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	45.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001153779</u>
Locus ID:	241197


[View online »](#)

UniProt ID: G3UYR4

RefSeq Size: 3314

Cytogenetics: 1 E2.1

RefSeq ORF: 1191

Synonyms: 9830131G07; BB233602; Serpinb10-ps