

Product datasheet for TP506376

Klhdc2 (NM_027117) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse kelch domain containing 2 (Klhdc2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206376 protein sequence Red =Cloning site Green =Tags(s) MADGNEDARAEDLPGPAFENYEAMELACPAERSGHVAVSDGRHMFVWGGYKSNQVRGLYDFYLPREELWI YYMETGRWKKINTEGDVPPSMMSGCAVCVDRVLYLFGGHHSRGNTNKFYMLDSRSADRGLQWERIDCQGI PPSSKDKLGVWVYKNKLIFGGYGYLPEDKVLGTFFDETSEFWNSSHPRGWNDHVVHILDTETFAWSQPIT TGKAPSPRAAHACATVGNGKGFVGGRYRDARMNDLHYLNLDTWENELIPQGVCPVGRSWHSLTPVSSDH LFLFGGFTTEKQPLSDAWTYCISKNEWIQFNHPYVEKPRLWHTACASDEGEVIVFGGCANNLLVHHRAAH SNEVLIFSVQPKSLVRLSLEAVICFKEMLANWSCLPKHLLHSVNQRFGSNNTSGS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	46 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: [NP_081393](#)
Locus ID: 69554
UniProt ID: [Q4G5Y1](#)
RefSeq Size: 1946
Cytogenetics: 12 28.77 cM
RefSeq ORF: 1218
Synonyms: 2310022K15Rik; D12Ert522e; HCLP-1; LCP
Summary: Represses CREB3-mediated transcription by interfering with CREB3-DNA binding.
[UniProtKB/Swiss-Prot Function]