

Product datasheet for TP506392

Tspyl4 (NM_030203) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse TSPY-like 4 (Tspyl4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206392 protein sequence Red =Cloning site Green =Tags(s) MNGVEGNNELSLANTTTTPSHASEDLDLKQDQGLQEETDTVREMEAAGEAGADGGASPDSEHCGPELCFR V AENSCAAAARGLEDAPSPSKGGDAPSAPVAADDSSKNGCQLEGPHSPAKPKALEACGAVGLGSQQMPGP K KTKEMTTTKCAISVATGKEGEAGAAMQEKKGLQKEKKVAGGGKEETRPRAPKINCMDLSLEAIDQELSNVN AQADRAFLQLERKFGRMRLHMQRRSFIIQNIPGFWVTAFRNHPQLSPMISGQDEDMRYMINLEVEEL K QPRVGCKFKFIFQSNPYFRNEGLVKEYERRSSGRVSLSTPIRWHRGQEPQAHHRNREGNTIPSFNWF SDHSLLEFDRIAIEIIGELWSNPLQYYLMGDGPRRGVRVPPRQPVESPRSRFQSG TR TRPLE QKLISEEDLA ANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	44.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: NP_084479
Locus ID: 72480
UniProt ID: Q8VD63
RefSeq Size: 3916
Cytogenetics: 10 18.85 cM
RefSeq ORF: 1218
Synonyms: 2610102M01Rik; B230210I21Rik; D10Bwg0791e