

## Product datasheet for TP506532

### Stk32b (NM\_022416) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse serine/threonine kinase 32B (Stk32b), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206532 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MGGNHSKPPVFDENEVNFDFHFQILRAIGKGSFGKVCIVQKRDTKKMYAMKYMKNQKCVERDEVNRNFR  
ELQIMQGLEHPFLVNLWYSFQDEEDMFMVVDLLLGGDLRYHLQQNVHFTEGTVKLYICELALALEYLQRY  
HIIHRDIKPDNILLDEHGHVHITDFNIATVLKGSSEKASSMAGTKPYMAPEVFQVYVDGGPGYSYPVDWWS  
LGTAYELLRGWRPYEIH SATPIDEILNMFKVERVHYSSTWCEGMVSLKLLTKDPESRLSSLRDIQSM  
TYLADMNWD AVFEKALMPGFVFNKGR LNC DPTFEEMILESKPLHKKKKRLAKHRSRSTKDCSPLNGH  
LQQCLETVRKEFIIFNREKLRRQQGHNGQLSDLDGRIGSQTSSKLDGRNNNLTHTCPRGCSS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-MYC/DDK
Predicted MW:	47.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_071861</a>
Locus ID:	64293
UniProt ID:	<a href="#">Q9JIX8</a>



[View online »](#)

RefSeq Size: 3458  
Cytogenetics: 5 B3  
RefSeq ORF: 1245  
Synonyms: 2510009F08Rik; Stk32; STKG6; YANK2