

Product datasheet for **TP506560**

Ankrd10 (NM_133971) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ankyrin repeat domain 10 (Ankrd10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR206560 protein sequence Red =Cloning site Green =Tags(s)
	MSAGAGVAVEAGFSSEELSLRFPLHRACRDGDLVALCSLLPHTPRAHLAAEDSFYGWTPVHWAHFGKL ECLMQLVRAGASLNVSTTRYAQTPAHIAAFGGHPQCLVWLIQAGANINKPDCEGETPIHKAARSGSLECI TALVGSGAHTDLRNASGLTAADIAQTQGFQECTQFLLSLQNHQMSRFCHNGTLSGGHESILPTHVSLGT N RKRCLEDSSESLGVKKARTRVPSLDHAMPLANGEAEDDADRMHVDRECAAVSDMKNSSSVLNTLTNGSVV N GHLDFFCANQLNGMESRSHPCLTGSNGVSNGQPLSSGQASVSANGTEEPEKTMGINPEMCGSLHLNGS PS SCVASRPSWVGDIGESLHYGHYHGFQDTAESIPELSSVLEHTSCVRVEQRYDSAVLGAMQLHHGS
	TR TRPLE QKLISEEDLA ANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	44 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_598732</u>
Locus ID:	102334
UniProt ID:	<u>Q99LW0</u>
RefSeq Size:	2291
Cytogenetics:	8 A1.1
RefSeq ORF:	1245
Synonyms:	4833425P12Rik; AW549277