

Product datasheet for **TP507313**

Tmem143 (NM_144801) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse transmembrane protein 143 (Tmem143), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207313 protein sequence Red =Cloning site Green =Tags(s)
	<p>MTVARWLRLLWERGQAMLHVTWGSKVRLWLSLVPALLGTPRALSSLENRMGVYRKMWNPKPCDWAQQYRER FIPFSKEQLLRLLIQEFHSSPAERAALFAFSAHVDFCTLFHYHQLLARLQALYDPINPDRETLDQPSLTD PERLSSEKDVQLALRPLLAQANFSPLESDALAYALVVHHPQDEVQVTINLDQYIYQFWALGQRVGMMPH KSSVSGSKRGFFRKLPPVERRYFKRVVLAARTKGGHVLKSKFDTPLEGLEQLLPELKV RTPV LQRALLNL MLVVS GVMIFVNVGMVILSDLKMATSLLLLLFAAFMGVKASKVFGQRRSAQALELAHVLYRSTSNNSL LSALALRAQEEHIKEALLAHSFLARRPGGSQGKPEETSRWLQSEVESWLLAQSGCDVTFNGPRALAHLQA LTPSLGLFPLPELPQLDPMVLGTPEATQAALGSSYPSP</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	51.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_659050
Locus ID:	70209



[View online »](#)

UniProt ID: [Q8VD26](#), [G3X9F4](#)

RefSeq Size: 2184

Cytogenetics: 7 B3

RefSeq ORF: 1377

Synonyms: 2310076O21Rik; AI481604