

Product datasheet for TP507482

Tmem151a (NM_001001885) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse transmembrane protein 151A (Tmem151a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207482 protein sequence Red =Cloning site Green =Tags(s) MPEGEGGDCGEVPALVPDGEPLREEQRPLKQSLGGSLCRESHWKCLLLTLIHACGAWVAWCRLATVPRL VLGPEAALARGAGGPPPTYPASPCSDGYLYIPLAFVSLLYLLYLAECWHCHVRSCQAPRTDANTVLALIH RLQQAAPCVWWKATSYHYVRRTRQITRYRNGDAYTTTQVYHERADSRTARGEFDYSAHGVRDVSKELVGL ADHAATRLRFTKCFSGSAEAEASYLTQRARFFSANEGLEDDYLEAREGMHLKDQDFRESLMVFADPRSP WYARAWVFWLVSAATLSWPLRVVAAYGTAHVHYQVEKLFGASSPPPGAVPSGPPLSRVATVDFTELEWHI CSNRQLVPSYSEAVVMGASSGAYLRGCQRCRRSVSSNSLPPARPSGPRLPFSRSLGAGGRTPGVFR SLGGPLGRRGEDTEPLESPPCYEDALYFPVLIVHGDSGCRGDGQAL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	51.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001001885</u>
Locus ID:	381199


[View online »](#)

UniProt ID: Q6GQT5

RefSeq Size: 4007

Cytogenetics: 19 A

RefSeq ORF: 1404

Synonyms: AW125688; Gm961; Tmem151