

Product datasheet for TP507565

Gbp111 (NM_029868) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse GC-rich promoter binding protein 1-like 1 (Gbp111), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207565 protein sequence Red =Cloning site Green =Tags(s) MAQHDFVPAWLNFTSPQSSKSSTATFDKHGEHLRGEGRFGISRRRHNSSDGFFNNGPLRTTGDSWHQ PS LFRHDSVDSGVSKGAYAGTTGWHGSSRGHDGMSQRSGGSSTGNHRHWNGSFHSRKGCAFQEKTPTEIR EE KKEDKVEKLQFEEEDFPSLNPEAGKQNQPCRPVGTSGVWENPPSAKQPSKMLVIKKISKEDPAAAFSAA FTSSGSHHANGNKVSTMVPSVYKNLVPKPAPPPSKPNAWKANRTEHKPGSLCSSRESAFTNPISVTKPVG LAAGAGLHSPKESPSSTTPIEISSRLTKLTRTTDRKSEFLKTLKDERNGDCSESRCDCKLEGLRLEG SHTPEPKENGEQGCLQNGLSLPMVEEREVLSSHLEAEHRLKAMGWQYEPENDESCLPLTEDELKEFHLR TEQLRRNGFGKNGLLQSRSCNLSPWRSTCIAECEDSDTETSSSQTSDDDAWK TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	52.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: NP_084144
Locus ID: 77110
UniProt ID: Q6NZP2
RefSeq Size: 3550
Cytogenetics: 4 D1
RefSeq ORF: 1419
Synonyms: 5330440M15Rik; BC002292
Summary: Possible transcription factor.[UniProtKB/Swiss-Prot Function]