

## Product datasheet for **TP507815**

### Suox (NM\_173733) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse sulfite oxidase (Suox), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207815 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MLLQLYRSVWVRLPQAIRVKSTPLRLCIQACSTNDSLEPQHPSLTFSDDNSRTRRWKVMGTLGLGWLV  
YHEHRCRASQESPRMYSKEDVRSNNPKTGVWVTLGSEVFDVTKFVDLHPGGPSKMLAAGGPPEPFWAL  
YAVHNQPHVRELLAEYKIGELNPEDSMSPSVEASDPYADDPHHPALRINSQRPFNAEPPPELLTEGYIT  
PNPIFFTRNHLPPVNLDPHTYRLHVVGAPGGQSLSLDDLHKFPKHEVTVTLCAGNRRSEMSKVKEVK  
GLEWRTGAISTARWAGARLCDVLAQAGHRLCDSEAHVCFEGLDSDPTGTAYGASIPARAMDPEAEVLLA  
YEMNGQPLPRDHGFPVRRVWVPGVVGARHVKWLGRVSVSEESYSHWQRDRDYKGFSPVDWDTVNFDLAPS  
IQELPIQSAITQPQDGAIVESGEVTIKGYAWSGGGRAVIRVDVSVDDGGLTWQEALEEGEEQCPRKAWAWR  
IWQLKAQVPAEQKELNIICKAVDSDSYNVQPDTVAPIWNLRGVLSNAWHRVHVQVVP

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-MYC/DDK
Predicted MW:	60.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_776094</a>



[View online »](#)

Locus ID: 211389

UniProt ID: [Q8R086](#)

RefSeq Size: 2349

Cytogenetics: 10 D3

RefSeq ORF: 1641

Synonyms: SO