

Product datasheet for **TP507938**

Pnkp (NM_021549) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse polynucleotide kinase 3'-phosphatase (Pnkp), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207938 protein sequence Red =Cloning site Green =Tags(s)

MSQLGSRGRLWLQSPGGPPPIFLPSDGQALVLRGRLTQVTRKCSRNVQVELIADPESRTVAVKQLGTVN
PSTVGVHELKPGSLGSLGDLVLYLVNGLYPLTLRWEELSTSGSQPDAPPDTPGDPEEGEDTEPQKKRVR
KSSLGWESLKKLLVFTASGVKPKQGVAAFDLDTLITRSGKVFPTSPSDWRILYPEIPKKLQELAAEGY
KLVFTNQMGIGRGKLPAEVFKGKVEAVLEKLGVPFQVLVATHAGLNRKPVSGMWDHLQEANEGIPISV
EDSVFVGDAAGRLANWAPGRKKKDFSCADRLFALNVGLPFATPEEFFLKWPAARFELPAFDPRTISSAGP
LYLPESSELLSPNPEVVAVGFPAGKSTFIQEHLVSAGYVHVNRTLGSWQRCVSSCAALRQGKRWVI
DNTNPDVPSRARIQCAKDAGVPCRCFNFCATIEQARHNNRFREMTDP SHAPVSDMVMFSYRKQFEPPTL
AEGFLEILEIPFRLQEHLDPALQRLYRQFSEG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	57.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_067524</u>



[View online >](#)

Locus ID: 59047

UniProt ID: [G5E8N7](#)

RefSeq Size: 2134

Cytogenetics: 7 B3

RefSeq ORF: 1569

Synonyms: 1810009G08Rik; PNK