

Product datasheet for TP507972

Arm9 (BC028964) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse armadillo repeat containing 9 (cDNA clone MGC:39015 IMAGE:5364429), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR207972 protein sequence Red =Cloning site Green =Tags(s) MADYHNLIGVTAELVDSLEATVSGKMITPEYLQSVCVRLFNSQMRQSLAHSVDFTPRGTASTMLRASLAP EKLKDVPLPLSLDYELKKDLIWGSDRLKAFLLQALRWRLTTSHPGEQRETVLQAYISNDLLDCHSHNQR SVLQLLSKSEAVRQYMARLINALASLAEGRLYLQNTKVLRLMLEGRLKEEDKDVITRENVLGALQKFSL RRPLQTAMIRDGLIFWLIDLLKDPDCLSDYTLSESVALLMNLCLRSAGKNMCAKVAGLMLKVLSDLLGHE NHEIQPYVNGALYSILSIPSIREEARAMGMEDILRCFIKEGNAEMIRQIEFIKQLNSEDLLDGVLESDD DEDEDDEEDHDIMEADLDKDELIQPQLGELSGEKLTTTEYLGIMTNTGKARRKGPASVQWSGDEPLRRPV TPGGHRTGCPVLGDHLISPNQAQQARNGCLRAMPVAHPDDYKEGKPGVTGCGTSSSCKNAVGAKPVLS SW AQESKR TRTRPLE QKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	55.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

Locus ID: 78795
UniProt ID: [Q9D2I5](#)
RefSeq Size: 1735
Cytogenetics: 1 C5
RefSeq ORF: 1488
Synonyms: 3830422A13Rik; 4831423D23Rik; 4930438O05Rik; 5730415N24Rik
Summary: Required for ciliogenesis.[UniProtKB/Swiss-Prot Function]