



**Locus ID:** 98053

**UniProt ID:** [Q3THK3](#)

**RefSeq Size:** 1720

**Cytogenetics:** 17 D

**RefSeq ORF:** 1527

**Synonyms:** 2810405L04Rik; C76800; RAP74

**Summary:** TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation (By similarity).[UniProtKB/Swiss-Prot Function]