

Product datasheet for TP508218

Depdc7 (NM_144804) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse DEP domain containing 7 (Depdc7), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208218 protein sequence Red =Cloning site Green =Tags(s) MATVREKAAALNLSALHSPAQRPPGFSVAQKPFGATYVWSSIINTLQTQVEVKRRHHLKRHND CFVGSE AVDVIYSHLTQNKFFGDVDIPRAKVVRVCQALMDYKVFEAVQTRVFGDKDKPTFEDSSCSLYRFTTIPNQ DNQLGHENKVHSPSRYSASFMSDDIKSDSLEDLWENLSLKPASSPHIHISTRLSPQVINEVWQEETIGR LLQLVDLPFLDSLLKQKEVVPKAPQSKRQPD MVNTSNYLDRGILRAYGDSQEDEWISAAIDCLEYLPDQM VVDISRNFP EQPDRTDLVKELLFDAIGRYSTREPLLNHLYDVHNGIAELLVNGKTEMALEATQLFLKLQ DSQNREEFRRLLYFMAVAADPSEFKLQEESDNRMVVKRVFSKAIVNNKNVSKGKTDLLVFLMDHQKDV KIPAA LHKIVSVKLLAIQKGRDPN KDTGYIYCQRID ESEYSSSTKKKTKDQLLSLLKTI DEDSKLSAKEK KKLLGQFRKTHPDIFTEYFGD TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	58.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_659053</u>


[View online »](#)

Locus ID:	211896
UniProt ID:	<u>Q91WS7</u>
RefSeq Size:	1734
Cytogenetics:	2 E2
RefSeq ORF:	1533
Synonyms:	AV216087