

Product datasheet for **TP508491**

Nelfa (NM_011914) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse negative elongation factor complex member A, Whsc2 (Nelfa), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208491 protein sequence Red =Cloning site Green =Tags(s)

MASMRESDTGLWLHNKLGATDELWAPPSIASLLTAAVIDNIRLCFHRLSSAVKLLKLLGLTLHLPRRTVDE
MKAALMDIIQLATLSDPWVLMVADILKSFPDGTSLNLDLEEQNPVQDILGELREKVSECEASAMLPLE
CQYLNKNALTTLAGPLTPPVKHFQLKRKPKSATLRAELLQKSTETAQQLKRSAGVPFHAKGRGLLRKMDT
TTPLKGIPKQAPFRSPTTPSVFSPSGNRTPIPPSRTPLQKERGVKLLDISELNTVGAGREAKRRRKTLDT
EWVEKPTKEETVENATPDYAAGLVSTQKLGSLNSEPTLPSTSYLPSTPSVVPASSYIPSETPPAPPSR
EASRPPEEPSAPSPTLPTQFKQRAPMYNSGLSPATPAPAAPTSPLTPTTPPAVTPTAQTTPVAMVAPQTQ
APAPVQQQPKKNLSLTREQMFAAQEMFKTANKVTRPEKALILGFMAGSRENPCPEQGDVIQIKLSEHTED
LPKADGGSTTMLVDTVFEMNYATGQWTRFKKYKPMTNVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	57.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_036044</u>



[View online »](#)

Locus ID: 24116

UniProt ID: [Q8BG30](#)

RefSeq Size: 2162

Cytogenetics: 5 B2

RefSeq ORF: 1593

Synonyms: Whsc2; Whsc2h

Summary: Essential component of the NELF complex, a complex that negatively regulates the elongation of transcription by RNA polymerase II (By similarity). The NELF complex, which acts via an association with the DSIF complex and causes transcriptional pausing, is counteracted by the P-TEFb kinase complex (By similarity).[UniProtKB/Swiss-Prot Function]