

Product datasheet for **TP508515**

Rcbtb1 (NM_027764) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse regulator of chromosome condensation (RCC1) and BTB (POZ) domain containing protein 1 (Rcbtb1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208515 protein sequence Red =Cloning site Green =Tags(s)
	MVDVGKWPIFTLLSPQEAGSIRKACVFGTSANEAIYVTDNDEVFVGLNYSNCLGTGDNQSTLVPKKLEA LCGKKIKLSYSGSPHVLLTTEDGVVYAWGHNGYSQLGNGTTNQGIAPVQVCTNLLIKQVIEVACGSHHS MALAADGELFAWGYNNCGQVGSSTANQPTPRKVTNCLHTKRNVNIACGQTSSMAVLDSGEVYGWGY NGN GQLGLGNNGNQLTTPRVAALHGMVCVNQIVCGYAHTLALTDEGLLYAWGANTYGQLGTGSKNNLLSPTQ IM VEKERVIEIAACHSTHTSAAKTQGGHVYMWGQCRGQSVILPHLTHFCCTDDVFACFGTPAVSWRLLSVEH EDFLTVAESLKKEFDSPEIADLKFRIDGKYIHVHKAVLKIRCEHFRSMFQSYWNEDMKEVIEIDQFSYPV YRAFLQYLYTDTVDLPPEDAIGLLDLATSYCENRLKRLCQHIIKRGITVENAFSLFSAAVRYDAEDLEEF CFKFCINHLTEVTQTAAFWQMDGPLLKEFIAKASKCGAFKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	58.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_082040</u>
Locus ID:	71330
UniProt ID:	<u>Q6NXM2</u>
RefSeq Size:	3949
Cytogenetics:	14 C3
RefSeq ORF:	1593
Synonyms:	5430409I18Rik; AW111883; CLLD7; CLLL7
Summary:	May be involved in cell cycle regulation by chromatin remodeling.[UniProtKB/Swiss-Prot Function]