

Product datasheet for TP508591

Cwf19II (NM_001081077) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse CWF19-like 1, cell cycle control (S. pombe) (Cwf19II), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208591 representing NM_001081077 Red =Cloning site Green =Tags(s)
	MAQKPLRLACGDVEGKFDVLFNRVRTIQKSGNFDLLLCVGNFFGSAQDAEWEEYKTGNKKAPIQTYVL GANNEETANYFQGADGCELAENITYLGRKGVFTGSSGLQIVYLSGTESLDEPVPAAHSFSPKDVSSLRTML CSASQFKGVDILLTSPWPKYVGSFGNSSGEVDTKNCGSALISSLAVSLKPRYHFAALEKSYIERLPYRNH VVLQESAQHATRFIALANVGNPEKKKYLYAFSIVPMKLMVAELVKQPPDV TENPYRDSGKQAAGGKHIP APQEEACQFFDLSEKQGRKRPSTGRDTRPPHAKQPRKPPQPPGPCWFCLASPEVEKHLVVNIGTHCYL ALAKGGLSDDHVLILPIGHYQSVVELSAEVVEVEKYKATLQRFFKSRGKRCVLFERNYRSHHLQLQVIP VPLSCCATDDIKDAFITQAQEQQIELLEIPEHSDIKQIAQPGAAYFYVELDTGEKLFHRIKKNFPLQFGR EVLASEAILNIPEKADWRQCQTSKDEEEALARRFRKDFEPDFTLDD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	60.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001074546</u>
Locus ID:	72502
UniProt ID:	<u>Q8CI33</u>
RefSeq Size:	3581
Cytogenetics:	19 C3
RefSeq ORF:	1611
Synonyms:	2610528C06Rik; A1854304; AV336991