

## Product datasheet for TP508636

### Cct5 (NM\_007637) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse chaperonin containing Tcp1, subunit 5 (epsilon) (Cct5), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208636 representing NM_007637 Red=Cloning site Green=Tags(s)

MASVGTLAFDEYGRPFLLIKDQDRKSRLMGLEALKSHIMAAKAVANTMRTSLGPNGLDKMMVDDKGDVTI  
TNDGATILSMMDVDHQIAKLMVELSKSQDDEIGDGTGWWLAGALLEEAQLLDRGIHPIRIADGYEQA  
ARIAIQHLDKISDKVLVDINNPELIQTAKTTLGSKVINSCHRQMAEIAVNAVLTVADMERRDVFELIK  
VEGKVGGRLEDTKLIKGVIVDKDFSHQPMPKKVDAKIAILTCFPEPPKPKTKHKLDVMSVEDYKALQKY  
EKEKFEEMIKQIKETGANLAICQWGFDEANHLLLQNGLPVAVRWVGGPEIELIAIATGGRIVPRFSELTS  
EKLGFAGVWQEISFGTTKDKMLVIEKCKNSRAVTIFIRGGNKMIIEEAKRSLHDALCVIRNLIRDNRVY  
GGGAAEISCALAVSQEADKPTLEQYAMRAFADALEVIPMALSENSGMNPIQTMTEVRARQVKESNPALG  
IDCLHKGSNDMQYQHVIETLIGKKQQISLATQMVRMILKIDDIRKPGESE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	59.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_031663</a>



[View online »](#)

<b>Locus ID:</b>	12465
<b>UniProt ID:</b>	<a href="#">P80316</a>
<b>RefSeq Size:</b>	1769
<b>Cytogenetics:</b>	15 B2
<b>RefSeq ORF:</b>	1623
<b>Synonyms:</b>	CCT-epsilon; Ccte; mKIAA0098; TCP-1-epsilon; TCPE
<b>Summary:</b>	Component of the chaperonin-containing T-complex (TRiC), a molecular chaperone complex that assists the folding of proteins upon ATP hydrolysis. The TRiC complex mediates the folding of WRAP53/TCAB1, thereby regulating telomere maintenance. As part of the TRiC complex may play a role in the assembly of BBSome, a complex involved in ciliogenesis regulating transports vesicles to the cilia. The TRiC complex plays a role in the folding of actin and tubulin.[UniProtKB/Swiss-Prot Function]