

Product datasheet for TP508705

Eif3d (NM_018749) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse eukaryotic translation initiation factor 3, subunit D (Eif3d), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208705 protein sequence Red =Cloning site Green =Tags(s) MAKFMTPVIQDNPSGWGPCAVPEQFRDMPYQPFSGKDRLGKVADWTGATYQDKRYTNKYSSQFGGGSQYA YFHEEDETFQLVDTARTQKTAYQRNRMFAQRNLRRDKDRRNMVQFNLQTLPKSAKQKERERIRLQKKF QKQFGVRQKWDQKSQKPRDSSVEVRSDWEVKEEMDFPQLMKMRYLEVSEPQDIECCGALEYDKAFFDRIT TRSEKPLRSIKRIFHTVTTTDDPVIRKLAKTQGNVFATDAILATLMSCTRSVYSWDIVVQRVGSKLFFDK RDNSDFDLLTVSETANEPQDEGENSFNSPRNLAMEATYINHNFSQQCLRMGRERYNFPNPNPFVEDDMDK NEIASVAYRYRRWKLGDIDLIVRCEHDGVMTGANGEVFSFINIKTLNEWDSRHCNGVDWRQKLDSSQRGAV IATELKNNYSKLARWTCCALLAGSEYLKLGYSRYHVKDSSRHVILGTQQFKPNEFASQINLSVENAWGI LRCVIDICMKLEEGKYLKDPNKQVIRVYSLPDGTFSSSEEEDEEEEEEEEEEEEEEEEE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	64 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_061219</u>



[View online »](#)

Locus ID: 55944

UniProt ID: [O70194](#)

RefSeq Size: 1910

Cytogenetics: 15 E1

RefSeq ORF: 1647

Synonyms: 66/67kDa; AA407891; eIF3p66; Eif3s7

Summary: mRNA cap-binding component of the eukaryotic translation initiation factor 3 (eIF-3) complex, a complex required for several steps in the initiation of protein synthesis of a specialized repertoire of mRNAs. The eIF-3 complex associates with the 40S ribosome and facilitates the recruitment of eIF-1, eIF-1A, eIF-2:GTP:methionyl-tRNAⁱ and eIF-5 to form the 43S pre-initiation complex (43S PIC). The eIF-3 complex stimulates mRNA recruitment to the 43S PIC and scanning of the mRNA for AUG recognition. The eIF-3 complex is also required for disassembly and recycling of post-termination ribosomal complexes and subsequently prevents premature joining of the 40S and 60S ribosomal subunits prior to initiation. The eIF-3 complex specifically targets and initiates translation of a subset of mRNAs involved in cell proliferation, including cell cycling, differentiation and apoptosis, and uses different modes of RNA stem-loop binding to exert either translational activation or repression. In the eIF-3 complex, EIF3D specifically recognizes and binds the 7-methylguanosine cap of a subset of mRNAs.[UniProtKB/Swiss-Prot Function]