

Product datasheet for TP508855

Msantd2 (NM_146222) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse Myb/SANT-like DNA-binding domain containing 2 (Msantd2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR208855 protein sequence Red =Cloning site Green =Tags(s)
	MAAPCGSELPANSPLKIPKMEVLSPASPGDLSGDNPSLSDPSTPRGASPLGPGSAAGSGAAASGGLGLGL GGRGAASSSVSFSPGGSGGAAAAAAACRGMSWTPAETNALIAVWGNERLVEARYQQLEGAGTVFGSKA PGPAMYERVSRALAEGLYERTPSQCRERIKTLRRCYSRVKEHGVGKRKSSYTFEQLQVFGQGWDAPC QPVLINSSGLYQELEDGSTMEDYSQEDWGNHSQELHGYPTDQELDEMPVSKRTLKIKQESSEEAQKRDT MQNIVQILESVQLKWELFQSWTDFSRHLHLSNKLAIFGIGYNTRWKEDIRYHYAEISSQVPLGKRLREYFN SEKPEGRIIMTRVQKMNWKNVYYKFLEITISEARCLEHMEIDWIPAHSKPTGGNVVQYLLPGGIPKSP GLYAIGYEECIERPLSPDVERHALDPGKEGRVDLETLSAQASLQVEVEPTRIIYCYLGIAEVRTLQQCLF LHFQANAKTFSKEWVGINGFLSQNCIVDPGVSPKSIYIKFVEVERDLSAGSLVECLEKAIGYPLKFNN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	61.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_666334</u>
Locus ID:	235184
UniProt ID:	<u>Q6NZR2</u>
RefSeq Size:	2307
Cytogenetics:	9 A4
RefSeq ORF:	1677
Synonyms:	2810450G17Rik; 9530092B10Rik