

Product datasheet for **TP509275**

Utp6 (NM_144826) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse UTP6 small subunit processome component (Utp6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209275 protein sequence Red =Cloning site Green =Tags(s)
	MAEIIQERIEDRIPELEQLERIGLFSHAEIKAIKKASDLEYKIHRRRTLLKEDFINYVQYEINLLELIQR RRARIKYSFKKDEIEYSMVHRVQGVFGRASAKWKDDVQLWLSYIVFCKKWGKTHTLSKIFSAMLAIHSNK PALWIMAAKWEMEDRLSSESARQLFLRALRFHPECPKLYQEYFRMELMHAELRKEKQEFKAAMDMDG DF DHPEEILKGELARIYKNSISKIKGAEFHVSLIAIQLFDFAKDLQKEIYDDLQALHTDDPLTWDYVARR ELEIESQPGEEQPVSKQAKAVEMGRREERCCAVYEEAVKALPTEAMWKCITFCLEFSSKTSVPLRGQ RLERTMLAFRKAHELKLLSEVQYKQWIDLLLRLQDLFKEALQVAEAGTELFKDSVTMWQTKLQVLIDSKSP DVEMRFEEAFAHLKPQVCLPLWISWAWESESAKSQEDTEAIFKKAIIVTGASSVTLKEKYLDWAYRSGG YKKARAVFKSLQESRPFSEFFRMMMQFEKEQEPCKMVNLREYYERALREFGTSDSDLWMDYIKEELNHP FGKPENCGQIYWRAMKMLQGQSAELFVAKHAMHQAGH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	70.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_659075</u>
Locus ID:	216987
UniProt ID:	<u>Q8VCY6</u>
RefSeq Size:	3871
Cytogenetics:	11 47.29 cM
RefSeq ORF:	1791
Synonyms:	4732497O03Rik; HCA66
Summary:	Involved in nucleolar processing of pre-18S ribosomal RNA.[UniProtKB/Swiss-Prot Function]