

Product datasheet for **TP509277**

Mettl25 (NM_207522) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse methyltransferase like 25 (Mettl25), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209277 protein sequence Red =Cloning site Green =Tags(s)
	MAAACPLPRTPDLPTLHDKLQGLLSFLRGALAISSAHTVDFYTKSVWQELVDLPPEVLAALRESAVEAE PREAETGSGFTELPKIFCETSQKLLSVEAFARTAKHYSVQNLGLCTPSEQLLTALQGNKRQRVDENVKAI EFMNTKKSHEVQAMSELICSIADYCGLKQIIDVGSGKGYLSSFLSLKYGLNLYGIDSSNTNTHGAKERNR KLKKHWSLYHPHSRADANGWASERPRELKVPKGVECKGDAESVQSRSLGNPDLSATDGLPDFSGSAISVI RKQQKNVLAQPAEEENLYFEDAFSLIDFLPVDAIEPTSSQVQNTKSGLRKERRNTASKARDSSIYSPLT SFITADSQLHDIEDLEDCLMVGLHTCGDLAPSTLRIFTSKAEVKAVCSVGCCYHLLSEEFENQHKDRCA NENWGPFMCHYLKEERWCCGRNARMSACLALQRVAVGGQLPTESLFYRAVLQNIKDYYGISKCEQHVG K IYKCSSFLEYVRMSLKKLGLDESKVSEEIIMDYENYKPRMNELEAFNMLKVVLPACIETLILLDRLCY LKEQDGVAWSALVKLFDPVQSPRCYAVIALKKQCDLG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	66.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_997405</u>
Locus ID:	216292
UniProt ID:	<u>Q6NXH8</u>
RefSeq Size:	2039
Cytogenetics:	10 D1
RefSeq ORF:	1791
Summary:	Putative methyltransferase.[UniProtKB/Swiss-Prot Function]