

Product datasheet for TP509301

Eps8l3 (NM_133867) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse EPS8-like 3 (Eps8l3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209301 protein sequence Red =Cloning site Green =Tags(s) MSRPSSRAIYLHRKEYSQSMASEPTLLQHRVEHLMTCCLGTQRVREPKDALQKLQEMDAQGRVWSQDL FL QVRDGLHLLDIETKEELDSYRLDNIKAIDVALNTCSYNSILSVTVQESGLPGISTLLFQCQEVGAELR TSLQKALEEELEERPRFGVHHPSQDRWKGPPLEPLPIQQAPPLEQRFSPPEHRFPPEQPHNMTSERSISP SSRSLTHYPSAREPNGFTLPPPPRRAPSPEDPERDEEVLNHVLRDIELFAGKLKEVQARNSHKTKLGRK KKKSKNGITQAEYIDCFQKIKLSFNLLGKLALRMQETSAPEFVGLIFQTLKFILSQCEAGLPKVISPL LTPKAIDLLQSCLSPPEDTLWKS LGTSWTTSWADWTGSEPPPYQPTFYDGWQIPQPRSMMPITNQDSISL RGSRMRSSLHFP RDEPYNHNPEYEDSNLPLSSPSPGRAALKMQVLYEFEARNAQELTVAQGEILEVLDQS KRWWLVKNEAGLTGYIPSNILEPLPAGAPRGHRQPSFRAPMLRLSSKPEEVTAWLQAENFSTVTVRTLGS LMGSQLLHMRPGELQMLCPQEAPRIQARLDAVRRLGMTH TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	68.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq:	<u>NP_598628</u>
Locus ID:	99662
UniProt ID:	<u>Q91WL0</u>
RefSeq Size:	2313
Cytogenetics:	3 F2.3
RefSeq ORF:	1800
Synonyms:	AI504489; Eps8r3