

Product datasheet for **TP509409**

Btbd9 (NM_172618) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse BTB (POZ) domain containing 9 (Btbd9), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209409 protein sequence Red =Cloning site Green =Tags(s)
	MSNSHPLRPFTAVGEIDHVVHILSEHIGALLIGEEYGDVTFVVEKKHFPAAHRVILAARCQYFRALLYGGMR ESQPEAEIPLQDTTAEAFMMLLRYIYTGRTLTDEKEEVLLDFLSLAHKYGFPELEDSTSEYLCTILNIQ NVCMTFDVASLYSLPKLTCMCCMFMDRNAQEVLASDGFLSKTALLNIVLRDSFAAPEKDIFLALLNWC KHNAKENHAEIMQAVRLPLMSLTELLNVVRPSGLLSPDAILDAIKVRSESRMDLNYRGM LIPEENIATM KYGAQVWKELKSALLDGGDTQNYDLDHGFSRHPIDDDCRSGIEIKLGQPSIINHIRLLLWDRDSRSYSYF IEVSMDELWDWIRVIDHSHYLCRSWQKLYFPARVCYIRIVGTHNTVNKIFHIVAFECMFTNKAFTLEKGL IVPLENVATIADCASVIEGVSRSRNALLNGDTKNYDWDSGYTCHQLGSGAIVVQLAQPYIIGSIRLLLWD CDDRSYSYYVEVSTNQQQWTMVADRTKVSCSWQSVTFERQPASFIRIVGTHNTANEVFHCVHFECPEQ Q SNQKEDSSEEPGTGDPSTPNQQLDPHAPRAPSASSLPPSPGPNSRSPNQNNQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	69.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_766206</u>
Locus ID:	224671
UniProt ID:	<u>Q8C726</u>
RefSeq Size:	6657
Cytogenetics:	17 A3.3
RefSeq ORF:	1836
Synonyms:	1700023F20Rik; 4930402L05; AI464421