

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_663331
Locus ID:	207259
UniProt ID:	Q8VCZ7
RefSeq Size:	4591
Cytogenetics:	18 E3
RefSeq ORF:	1857
Synonyms:	B230208J24Rik; Kr-pok; Zbtb36
Summary:	May be a tumor suppressor gene.[UniProtKB/Swiss-Prot Function]