

Product datasheet for **TP509507**

Tmem184c (NM_145599) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse transmembrane protein 184C (Tmem184c), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209507 representing NM_145599 Red =Cloning site Green =Tags(s)
	MPCACNRSNWRRWIRPLLVLFYATTILVAVPICIWKFKMKVGMHTKSWFIAGIFLLLTIPVSLWGILQH LVHYTQPELQKPIIRILWMVPIYSVDSWVALVYPKIAIYVDTWRECYEAYVIYNFMIFLTNYLTIRFPNL ILHLEAKDQQNHILPLCCPPWAMGEMLLFRCKLGVLYQYTVVRPITVTALVCEILDVYDEGNFGFSNAW TYLVILNNLSQLFAMYCLLLFYKVLKEELSPIQPVGKFLCVKLVVVFSFWQAVLIALLVKLGVISEKRTW EWQSAEAVATGLQDFIICIEMFFAAIAHHYTFSYKPYVHEAEEGSCFDSFLAMWDVSDIRDDISEQVRRV GRTMRGYPKKKCFPGDPDHNEHSSLLSSSSQDLTSGSSKVPSPGGLYQGFGHTISSQSPISIASIYEEIM NDIPEEQKLLNPGQDVTINIPPEEQKLNLTGKDVMLDRQTITSPRPISIASIYEKIMIDIPEEQKLTDT GQDVMINRQTITPPRPISITSIYEEVLIDITEKQKLNLTGKNVMINIPEEQKLNIDKRKDVMIPEEQN AIPDNSQYQDQEQIVTLQALFPSTETSENSMIDTSESQEQSSDLCTESSDSSTESSDLSTDP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	71.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_663574
Locus ID: 234463
UniProt ID: Q3TPR7
RefSeq Size: 2701
Cytogenetics: 8 C1
RefSeq ORF: 1866
Synonyms: 8430433H16Rik; AI645453; AL033298; BC004056; Tmem34
Summary: Possible tumor suppressor which may play a role in cell growth.[UniProtKB/Swiss-Prot Function]