

## Product datasheet for TP509617

### Hnrnpr (NM\_028871) Mouse Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse heterogeneous nuclear ribonucleoprotein R (Hnrnpr), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA</b>	>MR209617 representing NM_028871
<b>Clone or AA</b>	Red=Cloning site Green=Tags(s)
<b>Sequence:</b>	<p>MANQVNGNAVQLKEEEPMDTSSVTHTEHYKTLIEAGLPQKVAERLDEIFQTGLVAYVDLDERAIDALRE FNEEGALSVLQQFKESDLSHVQNKSAFLCGVMKTYRQREKQGSKVQESTKGPDEAKIKALLERTGYTLDV TTGQRKYGGPPPDSVYSGVQPGIGTEVFVGKIPRDLYEDELVPLFEKAGPIWDLRLMMDPLSGQNRGYAF ITFCGKEAAQEAVKLCDSYEIRPGKHLGVCISVANNRLFVGSIPKNKTENILEEFSKVTEGLVDVILYH QPDDKKKNRGFCFLEYEDHKSAQAARRRLMSGKVKVWGNVTVWADPVEEPDPEVMAKVKLVFVRNLAT TVTEEILEKSFSEFGKLERVKKLDYAFVHFEDRGA AVKAMDEMNGKEIEGEEIEIVLAKPPDKKRKERQ AARQASRSTAYEDYYYHPPPRMPPPMRGRGRGGRGGYGYPPDYGYEDYDDYGYDYHDYRGGYEDPYY GYDDGYAVRGRGGGRGGGAPPPRGRGAPPPRGRAGYSQRGAPLGPGRGSRGGRGGPAQQQRGRGSRGA RGNRGGNVGGKRKADGYNQPDSKRRQTNNQNWGSQPIAQQPLQQGGDYSYGYNNNDNQEFYQDTYGGQ WK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	71.3 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_083147](#)  
Locus ID: 74326  
UniProt ID: [Q8VHM5](#)  
RefSeq Size: 8025  
Cytogenetics: 4 D3  
RefSeq ORF: 1896  
Synonyms: 2610003J05Rik; 2610528B01Rik; Hnrpr