

Product datasheet for TP509671

Cttnbp2nl (NM_001163332) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse CTTNBP2 N-terminal like (Cttnbp2nl), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209671 protein sequence Red =Cloning site Green =Tags(s) MNLEKLSKPELLTLFSILEGELEARDLVIEALKAQHRDTFIEERYGKYNISDPLMALQRDFETLKEKNDSEKQPVCTNPLSVLKAVMKQCKNMQERMLSQLAAAESRHRKVILDLEEERQRHAQDTAEGDDVTYMLEKE R ERLTQQLEFEKSQVKKFEKEQKKLSSQLEEERTRHKQLSSMLVLECRKATSKAAEEGQKAGELSLKLDKEKSRASKLEEELAAERKRGLQTEAQVEKQLSEFDIEREQLRAKLNREENRTRALKEEVESLKKLVKDLEAAQQHRSTSEQGREPVTMSRGTEPPMRVSAFCQTESVQTERSHGSVITKLTDTGLPGPTTAAYSYAKANGHCDPEIQTRELTSDDSSTENQGPPREKSAAVAQEKPVENGGCPVGTETPVTMPSHLPSSGSSLSPSSTASSLTSSPCSSPVLTKRLLGSAASSPGYQSSYQVGINQRFHAARHKFQSQADQDQQASGLQSPPSRDLSPDLLDNSAAKQLARNTVTQVLSRFTNQGPVKPVSPNSSPFGTDYRNLASTASPRGDTSHSPTPGKVSSPLSPLSPGIKSPTIPRAERGNPPPIPPKKPGLTPSQSATTPTVTKTHSQASSLAATEDLASSCSPSAVVANGKDV EILLPTSS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	69.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001156804</u>
Locus ID:	80281
UniProt ID:	<u>Q99LJ0</u>
RefSeq Size:	4838
Cytogenetics:	3 F2.2
RefSeq ORF:	1914
Synonyms:	AA552995; AA589392; AU018624; BC003236; mKIAA1433
Summary:	Regulates lamellipodial actin dynamics in a CTTN-dependent manner.[UniProtKB/Swiss-Prot Function]