

Product datasheet for **TP509764**

Sparcl1 (NM_010097) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse SPARC-like 1 (Sparcl1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209764 protein sequence Red =Cloning site Green =Tags(s)
	<p>MKAVLLLLCALGTAVAIPTSTRFLSDHSNPTTATLVTPEDATVPIAGVEATADIENHPNDKAEKPSALNS EEETHEQSTEQDKTYSFEVDLKDDEEDGDGDLSDVPTGTLTLDLQEGTSEPQQKSLPENGDFPATVSTSY VDPNQRANITKGKESQEQPVSDSHQQPNESSKQTQDLKAEESQTQDPDIPNEEEEEEEEEEEEEEEEEPEP IGAPSDNQEEGKEPLEEQPTSKWEGNREQSDDTLEESSQPTQISKTEKHQSEQGNQGGQESDSEAEGEDKA AGSKEHIPHTEQQDQEGKAGLEAIGNQKDTDEKAVSTEPTDAAVVPRSHGGAGDNGGGDDSKHGAGDDYF IPSQEFLEAERMHLSYLYKYGGEETTTGESENRRRAADNQEAKKAESSPNAEPSDEGNSREHSAGSCT NFQCKRGHICKTDPQGKPHCVCPETCPPAKILDQACGTDNQTASSCHLFATKCRLEGTKKGHLQLD YFGACKSIPACTDFEVAQFPLMRDWLKNILMQLYEPNPKHGGYLNEKQRSKVKKIYLDEKRLLAGDHPI ELLLRDFKKNYHMYVPVHWQFNELDQHPADRILTHSELAPLRASLVPMEHCITRFFEECDPNKDKHITL KEWGHCFGIKEEDIDENLLF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	72.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_034227](#)
Locus ID: 13602
UniProt ID: [P70663](#), [Q3TYD5](#), [Q3UYK7](#)
RefSeq Size: 3132
Cytogenetics: 5 50.55 cM
RefSeq ORF: 1953
Synonyms: Ecm2; hevin; mast9; Sc1