

Product datasheet for TP509921

Lrrc45 (NM_153545) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse leucine rich repeat containing 45 (Lrrc45), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209921 protein sequence Red =Cloning site Green =Tags(s)
	MEEFRRSYNRLCEESGAEPQEAVLQQLHQLPKGGDLTTQSLTVETCRALGKLLHKETLLKELVLSDCML SEEGSTLLFQGLCANTSVQHLDLKGNNLRATGAEALGKLLRQNKSIQSLTLEWNNLGTWEDAFATFCGGL AANSALRQLDLRNNQISHKGAEELALALKGNTTLQQLDLRWNNIGLLGGRALVNCPLSNRTLWKLDLAG N NIPGDILRAVEQAMDHNQDRLTAFRENQARTKILSKEVQHLQEEKSKQFLDLMETIDKQREEMARDSRAS AVRVGQLQEALNERQSIINALKAKLQMTEAALALSEQKVRDLGELLVAGDQERQSLSQRHEKERKLERQE AADRESKLLRDLAANEKNLLRSQVDELERKARSQQEQFLTKQELTNTSAELKIRAIQAEERLDVEKR RAKQNMEDLEKLHSKEVDHMRHLEESERAMQERVQRLEALRLSLEEELSRMKA AVL SERGQAEELIKA RNQARLEEQHRLAHLEEKIRLLAQARDEAAGTCVQQKQMVAESQARVSQLNLQMEGQQRRLEELQQEL IN KDQEKVAEVARVRVELQEQMGRMQADLVAQEALREKVAALERQMKVIGSEHREALLDRESENASLREKLR LKEAEISRIRDEEAQRASFLQNAVLAYVQGSPLRALSPPK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	76.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



View online »

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_705773</u>
Locus ID:	217366
UniProt ID:	<u>Q8CIM1</u>
RefSeq Size:	2476
Cytogenetics:	11 E2
RefSeq ORF:	2010
Synonyms:	BC023296
Summary:	Component of the proteinaceous fiber-like linker between two centrioles, required for centrosome cohesion.[UniProtKB/Swiss-Prot Function]