

Product datasheet for **TP509986**

Dhx35 (NM_145742) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse DEAH (Asp-Glu-Ala-His) box polypeptide 35 (Dhx35), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209986 protein sequence Red =Cloning site Green =Tags(s)
	MAAPLGPVKFWKPGTEGPGVSISEERQSVTENS GTTIIYNPYAALSIEQQRQKLPVFKLRNHILYLVENY QTVVIVGETGCGKSTQIPQYLAEAGWTAEGRVVGVTQPRRVAAVTVAGRVADERGAVLGHEVG YCIRFDD CTDPLATRIKFLTDGMLVREMMVDPLLTKYSVIMLDEAHERTLYTDIAIGLLKKIQKKRGDLRLIVASAT LDAEKFRDFFNQNETSDPARDTSVTLTVEGRTFPVDIFYLQSPVPDYIKATVDTVVKIHQTEGDGDILAF LTGQEEVETVSMLEIQAARALARTGMKKHLRVLP MYAGLPSFEQMKVFERVSHSVRKVIVATNVAETSIT ISGIVYVIDCGFMKL RAYNPRTAIECLVVVPVSQASANQRAGRGRNRSGKCYRLYTEEAFDQLPQSTVP EMQRSNLAPVILQLKALGIDNVLRFHFMSPPPAQSMVQALELLYALGGLDKDCRLTEPLGMRIA EFPLNP MFAKMLLESGNFGCSQEVL SIAAMMQIQNVFVFP SNQKSQATRVHRKF AVEEGDHLTMLNVYEAFIKHN K SSQWCQEHLN YKGLVRAATVRDQLKKLLVKFQVPKVSSEGD PDPVLR C IVSGFFANAARFHSSGAYRVI YNEVIQTSKY YMRDVTAIESAWLLELAPHFYQQGTHLSLKAKRAKVQDQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	75.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_665685
Locus ID:	71715
UniProt ID:	Q8K1G9
RefSeq Size:	3218
Cytogenetics:	2 H1
RefSeq ORF:	2037
Synonyms:	1200009D07Rik; Ddx35