

Product datasheet for TP509994

Ushbp1 (NM_181418) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse USH1 protein network component harmonin binding protein 1 (Ushbp1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR209994 protein sequence Red =Cloning site Green =Tags(s) MSARATRPRSRGRHPLPGELDPVAESSEEDANNGCFETKPELHQEYLGPGQLQSPRNWIKEGCDVGS A VLSLSPEEQCEPEVEAHQVLQEAPLDSGPAEISVPSVYETLQCRSSLEAVVAALRHHSLSFPKSVEAED RDQGAPGPGFDEKEDAGPGQQAARLIERNAWLRLALCNREDELACTQASLQDAQAEKETLQRQVQELE D SLMQMEASPPTPILRAGRNRNSNSTSGAERRPWVPQDSPMAHPFLQRLRSDSSTQSFGLSTQHPAPET Y LMEDQMGQLQGSIEKLKCFNRLLLAVLQGYKGRCESLSIKLAKREAEATALRLALQYSEDCEEVYEVLLA LKMAGLGTGVGTTKGDQAEEKEASRLLVKKEATMEVETPQPSSGSGVDKPTPEELAVQLHGYVQHLRE RWALVKIPEELGPVTAPKATMPHAEATVQAILEIQPGPTLPHLEKSQIQDLAATRDSLADLVRLQLTQ REKRGLLELREAAALRAQGPAHQLLLQQLRWERAHFAGDGSSGGSSEDPSSDEEAGEDRQQHYQGPLALLD G QMGKVWDSENVSQELSASLARAIDLRAQLQSLRRQLEQVTQKGRTRRVQSAELSRELCKAHSALALAFRG AHRKQEEQRRKLEQQVARLQAQQAEEELAVLTATARALGKPGAPQPAQTFL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	74.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_852083
Locus ID:	234395
UniProt ID:	Q8R370
RefSeq Size:	3549
Cytogenetics:	8 B3.3
RefSeq ORF:	2040
Synonyms:	2210404N08Rik; 9430064H24; AIEBP; MCC2