

Product datasheet for TP510005

Rabl6 (NM_001024616) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RAB, member RAS oncogene family-like 6 (Rabl6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR210005 protein sequence Red =Cloning site Green =Tags(s) MFSALKKLVGSEQAPGRDKNIPAGLQSMNQALQRRFAKGVQYNMKIVIRGDRNTGKTALWHRLQGKKF VE EYIPTQEIQVTSIHWNKYKTTDDWVKVEVDVVDKGKCKKRGDGLKTENDPQEAESLALDAEFLDVYKNC NGVMMFDITKQWTFNYVLRELKVPVHPVPCVLGNYRDMGEHRVILPDDVRDFIEHLDRPPGSSYFRYA ESSMKNSFGLKYLHKFFNIPFLQLQRETLLRQLETNQLDIDATLEELSVQQETEDQNYSIFLEMMEARSR GHASPLAANGQSPSSGSQSPVPPSAVSTGSSSPSTQPAPQLSLGVSSTCPSAPSPVPSLEAMPSSVHS SAPAPTAPAPAPAQRRIISRLFGTSPAAEVTPSPPEPAPALEAPARVQNVEDFVPEDGLDRSFLEDTS VPKDKKVGAKGPQQDSDDGEALGGNPMVAGFQDDVDIEDQTHGKSLLPSDPMPSPKSNISLSSEEEAEG L AGHPRVAPQQCSEPETKWSSTKVSHQPQKKRAPTRGTPPWSGLTTDDSERPQEGKDKQVSSESDEGPIA AQMLSFVMDDPDFESDESDTQRRMGRFPVREDLSVTDVDEDTGPAQPPPPSKLPGAFLKNDSDLFGLGL E EMGPKESSDEDRDLSKPSKEKKKKKKKSKEEEEKTTKKKSKHKKSKDKEEGKEDRKKKRKPPRSKEQKAA DELEAFLGGGAPGSRHPGGGDYEEL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	79.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001019787
Locus ID:	227624
UniProt ID:	Q5U3K5
RefSeq Size:	3173
Cytogenetics:	2 A3
RefSeq ORF:	2175
Synonyms:	B230208H17Rik; parf; Rbel1; Rbel1a; Rbel1b
Summary:	May enhance cellular proliferation. May reduce growth inhibitory activity of CDKN2A (By similarity).[UniProtKB/Swiss-Prot Function]