

## Product datasheet for TP510049

### Samd4b (NM\_175021) Mouse Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse sterile alpha motif domain containing 4B (Samd4b), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>MR210049 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MMFRDQVGILASWFKGWNECEQTVALLSLLKRVTRTQARFLQLCLEHSLADCNDIHLLESEANSAAIVSQ  
WQQESKEKVVSLLSHLPLLQPGNTEAKSEYMRLQKVLAYSIESNAFIEESRQLLSYALIHATTLEDRL  
NALALWLSHLEERLASGFRTRPEPSYHSRQGSDEWGGPAELAPGEAGPGWQDKPPRENGHVPHFPSSVSP  
PAINSIGSNANTGLPCQIHPSPLKRSMALIPTSPQAPGEWPSPEELGARAAFTTPDHAPLSPQSSVASSG  
SEQTEEQGSSRNTEFQEDGSGMKDVPSWLKSLRLHKYAALFSQMSYEEMMTLREQHLESQNVTKGARHKIA  
LSIQKLRERQSVLKSLEKDVLEGGNLWNLALQELQIIITPIKAYSVLQATPTAKDEGRGPELLPGAEPPL  
THPGSDKGTEAKDPPAAENYPPPPAPAPSDSSEPAPAPVADGDIPSQFTRVMGKVCTQLLVSRPDEENIT  
SYLQLIEKCLTHEAFTETQKKRLLSWKQVQLKLLRTFPRKAALDMQSYRQKQGWAFGNSLPIAGSVGMG  
VARRTQRQFPMPPRALPPGRMGLLSPSGIGGVSPRHALTSPSLGGQGRQNLWFANPGGSNSMPSQSRSSV  
QRTHSLPVHSSPQAILMFPPDCVPGPDLEINPTLESLSLCLSMTEHALGDGTDKTSTI

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	75 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_778186](#)  
Locus ID: 233033  
UniProt ID: [G5E8A7](#)  
RefSeq Size: 4391  
Cytogenetics: 7 B1  
RefSeq ORF: 2064  
Synonyms: BC042901