

Product datasheet for **TP510096**

Zswim3 (NM_178375) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse zinc finger SWIM-type containing 3 (Zswim3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR210096 protein sequence Red =Cloning site Green =Tags(s)

MELGSCFKTYEDFKECFSAykkETRCsfivRDCISvRFHNLNNGTSFREDILYVQVKFVCIRTQSNRKKTRKVDRCPAYLLLQYNETLDRLFISELNTQHIHADSNASGGIPASKPQAICLHKLPPVQLSIRKDLDTAEKPSVEPSFLDKIHKPSNPEQEGISPSDMAKIAKVMKNFLTVDDEGSMASFSVGTSDQLDRLSFQSSKMSDLFARFPELTLHHRVENSQGHILYAFLVENKERESRVVHFAVLKAETAPSVAKMLNIFTEFNSDWPVKVMVFVDPSFPHRAILQEIFPSARTLLSIYHTTRLLEKKLHQSSVNASFKRLMKEALREAVFVSSDASLKNLCESQVLLDEELFSLQAHWFSCCELLWYIHVRKGLHACNTYMDSLIVTSKVVSLFREQRSLDCILHFVDYIDFFNTKGLKSLPTNAPKLRRTRLPSTPPRPKPFRCICGGDTRLPVEEVEETKADSAQSQLPQPQDQSSKAGMLDILHQSGSELAYKLCHNEWEVQNSTHLVDEAGSSVAVQLENSHQVSKDGCSCSCSFQQCYHLPCRHILALLHTSQQPVGAMVCCRWQKKYQHLLDPNGELQDRGIIPNTDQPEKQGQNHMIRDLSRELANLLMQTEGPELEERCSTLRKIVDIWADPYQLPESSQQPVDFRDVGCLPFLWGWKLEEDSPPLAEAMVHD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	79.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_848462](#)
Locus ID: 67538
UniProt ID: [Q8CFL8](#)
RefSeq Size: 2711
Cytogenetics: 2 H3
RefSeq ORF: 2088
Synonyms: 4921517A06Rik; C86566