

## Product datasheet for TP510114

### 2700050L05Rik (BC046235) Mouse Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse RIKEN cDNA 2700050L05 gene (cDNA clone MGC:58090 IMAGE:6817045), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>MR210114 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MADSDENGSYSTSSDPADDNKAVAIKSVGELSVPEKYKSIHQIRPSCAFPVCHDTEERCRLVLSYVLEG LKSVDSSIKKESDLPAADPSTPIPLKYEDSTRGGPEGLEKQMALFLDKMGSIQKGCSCGQSGMTPGSGWQ HKMKLQLILKSSKAYIILSDAAMSLQKYGRALRYIKLALQSHDITYCCLCTNMLSEVLLFLSQYLTLCGDI QLMLAQNANNRAAHLEEFNYQTKEDQEILHSLHRESSCQGFATDLSTDLESQSVSCKCYEAANEILQ FSDLKSQNPEHYVQVLKRMGNIRNEIGVFYMNQAAALQGERVVSKSVSAAEQQLWKKSFCFEKGIHNFE SIDDATNAALLCNTGRLMRVCAQAHCGAEDEFKREFSPEEGLYYSKAVDYLLKALRSLGTRDMHPVWD SVNWELSTTYFTMATLQQDYAPLSRKAQEIEKEVSEAMMKSLLKHCVDVSATARQPLCQYRAATIHHLRA SMYHSLRNQVGDEHLRQHRVLADLHYSKAAKLFQLLKDAPCELLRVQLERVAFAEFQMSSQNSNVGKL KTLGALDIMVRTEHAFQLIRKELVEECDQPKNDEATPAADSSPNLNREEVIKLLSIFESRSLFLLQSI KLMSSSKRKMSNAEEDIVLQTNKQIYSQLLRATANRNSTLLERIEVICLLEQLASGSSRSSGSAVP
	<b>TRTRPLEQLISEEDLAANDILDYKDDDDKV</b>
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	78.4 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>Locus ID:</b>	214764
<b>UniProt ID:</b>	<a href="#">Q6GQV7</a>
<b>RefSeq Size:</b>	2870
<b>Cytogenetics:</b>	7 F3
<b>RefSeq ORF:</b>	2094
<b>Synonyms:</b>	MGC27934, MGC58090
<b>Summary:</b>	Transcription factor involved in erythroid differentiation. Involved in transcriptional activation of the globin gene (By similarity).[UniProtKB/Swiss-Prot Function]