

## Product datasheet for TP510820

### Mcm6 (NM\_008567) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse minichromosome maintenance complex component 6 (Mcm6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR210820 protein sequence Red=Cloning site Green=Tags(s)

MDLAAAEPGAGSQHPEVRDEVAEKCKLFLDFLEEFQGS DGEIKYLQFAEELIRPERNTLWVSFADLEQ  
FNQQLSTTIQEEFYRVYPYLCRALKTFVKDRKEIPFAKDFYVAFQDLPTRHKIRELTSSRIGLLTRISGQ  
VVRTHPVHPELVSGTFLCLDCQTVIKDVEQQFKYTQPNICRNPVCANRKRFLDNTKSRFVDFQKVRIQE  
TQAELPRGSIPRSLEVILRAEAVESAQAGDRCDFTGALIVPVDVSKLSTPGARAETNSRVSGADGYETEG  
IRGLRALGVRDLSYRLVFLACHVAPTNPFRGGKELRDEEQTAESIKNQMTVKEWEKVFEMSQDKNLYHNL  
CTSLFPTIHGNDDEVKRGVLLMLFGGVPKTTGEGTSLRGDINVCIVGDPSTAKSQFLKHVDEFSPRAYTS  
GKASSAAGLTAAVRDEESHEFVIEAGALMLADNGVCCIDEFDKMDMRDQVAIHEAMEQQTISITKAGVK  
ATLNARTSILAAANPVSGHYDRSKSLKQNINLSAPIMSRFDLFFILVDECNEVTDYAIARRIVDLHSRIE  
ESIDRVYSLDDIRRYLLFARQFKPKISKESEDFIVEQYKRLRQRDGGSGVTKSSWRITVRQLESMIRLSES  
MARMHCCDEVQPKHVKEAFRLLNKSIIRVETPDVNLQEEEEIQMETDEGQGGVNGHADSPAPVNRFGSS  
EDASQETVSKPSLRLGFAEYCRISNLIVLHLRKMEEEEDESALKRSELVNWYLKEIESEIDSEELINKK  
TIIKVVHRLTHYDHVLIELTQAGLKGSSSEGSSESYEEDPYLVWNPNYLLED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	92.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online >](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_032593</a>
<b>Locus ID:</b>	17219
<b>UniProt ID:</b>	<a href="#">P97311</a> , <a href="#">Q542I2</a>
<b>RefSeq Size:</b>	2901
<b>Cytogenetics:</b>	1 55.8 cM
<b>RefSeq ORF:</b>	2466
<b>Synonyms:</b>	ASP-I1; D1Wsu22e; Mcmd6
<b>Summary:</b>	Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells. The active ATPase sites in the MCM2-7 ring are formed through the interaction surfaces of two neighboring subunits such that a critical structure of a conserved arginine finger motif is provided in trans relative to the ATP-binding site of the Walker A box of the adjacent subunit. The six ATPase active sites, however, are likely to contribute differentially to the complex helicase activity.[UniProtKB/Swiss-Prot Function]