

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	237339
UniProt ID:	Q8BLB7
RefSeq Size:	4280
Cytogenetics:	10 A4
RefSeq ORF:	2574
Synonyms:	AI481284; MBT-1
Summary:	Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility (By similarity). Required for normal maturation of myeloid progenitor cells.[UniProtKB/Swiss-Prot Function]